

Product Datasheet

C1QBP Rabbit mAb [KD Validated] (orb1950430)

Description	Binds to the globular heads" of C1Q thus inhibiting C1 activation.,,,,1,,1,1,1,0,1,1701055481,0,0,0,0,0,0,1701055481,0,-1,-1,0,0 324741,57,60074,ORP1 Rabbit mAb,0,10,ORP1-Rabbit-mAb,ORP1 Rabbit mAb,,0,50ul,249,249,0,0,0,58,ABWAY,CY10417,CY10417,0,0,0,Yes,Rabbit Monoclonal Antibody,Binds phospholipids; exhibits strong binding to phosphatidic acid and weak binding to phosphatidylinositol 3-phosphate.,,,,1,,1,1,1,0,1,1701055481,0,0,0,0,0,0,1701055481,0,-1,-1,0,0 324742,57,60075,MUTA Rabbit mAb,0,10,MUTA-Rabbit-mAb,MUTA Rabbit mAb,,0,50ul,249,249,0,0,0,58,ABWAY,CY10441,CY10441,0,0,0,Yes,Rabbit Monoclonal Antibody,"Involved in the degradation of several amino acids, odd-chain fatty acids and cholesterol via propionyl-CoA to the tricarboxylic acid cycle. MCM has different functions in other species.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human C1QBP
Form/Appearance	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Recombinant
Antibody Type	Recombinant Antibody

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Dilution Range WB 1:500~1:2000 IHC 1:50~1:200

Expiration Date 12 months from date of receipt.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of C1QBP expression in HeLa cell lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)