

Product Datasheet

General Malondialdehyde (MDA) ELISA Kit (orb1952272)

| | |
|----------------------------|---|
| Catalog Number | orb1952272 |
| Category | Assays and Kits |
| Description | General Malondialdehyde (MDA) ELISA Kit |
| Tested applications | ELISA |
| Application notes | Detect Range: 31.2-2000 pmol/mL. Sensitivity: 10.2 pmol/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 μ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Note | For research use only |
| Expiration Date | Please enquire. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)