

Product Datasheet

Rat Interleukin 22 (IL22) ELISA Kit (orb1952396)

| | |
|----------------------------|--|
| Catalog Number | orb1952396 |
| Category | Assays and Kits |
| Description | Rat Interleukin 22 (IL22) ELISA Kit |
| Reactivity | Rat |
| Tested applications | ELISA |
| Application notes | Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.17 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum,plasma,tissuehomogenates,celllysates,cellculturesupernatesorotherbiologicalfluids. |
| Note | For research use only |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)