

Product Datasheet

General Acetylcholine (ACH) ELISA Kit (orb1952492)

Catalog Number	orb1952492
Category	Assays and Kits
Description	General Acetylcholine (ACH) ELISA Kit
Reactivity	All
Tested applications	ELISA
Application notes	Detect Range: 15.6-1000 pg/mL. Sensitivity: 6.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum,plasma,tissuehomogenates,celllysates,cellculturesupernatesorotherbiolog icalfluids.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)