

Product Datasheet

Human Anti-Adalimumab antibody ELISA kit (orb1952887)

Catalog Number	orb1952887
Category	Assays and Kits
Description	Human Anti-Adalimumab antibody ELISA kit
Target	Anti-Adalimumab
Reactivity	Human
Assay Type	Sandwich
Range	qualitative
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Storage	2 - 8 °C for 12 months
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)