

Product Datasheet

Clusterin Protein, Human, Recombinant (sCLU, His) (orb1953091)

Description	Clusterin Protein, Human, Recombinant (sCLU, His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 51.51 kDa and the accession number is NP_001822.3.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	51.51 kDa (predicted); 33.5 kDa (reducing conditions)
Uniprot ID	P10909
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)