

Product Datasheet

OX40/TNFRSF4 Protein, Human, Recombinant (His & Avi), Biotinylated (orb1953093)

Description	OX40/TNFRSF4 Protein, Human, Recombinant (His & Avi), Biotinylated is expressed in HEK293 mammalian cells with His and Avi tag. The predicted molecular weight is 23.34 kDa and the accession number is NP_003318.1.
Storage	-20°C
Tag	C-His-Avi
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	23.34 kDa (predicted); 38.6 kDa (reducing conditions)
Uniprot ID	NP_003318.1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)