

## Product Datasheet

# VEGFR3/FLT4 Protein, Human, Recombinant (His & Avi), Biotinylated (orb1953331)

|                          |                                                                                                                                                                                                                     |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | VEGFR3/FLT4 Protein, Human, Recombinant (His & Avi), Biotinylated is expressed in HEK293 mammalian cells with His and Avi tag. The predicted molecular weight is 87.80 kDa and the accession number is NP_002011.2. |
| <b>Storage</b>           | -20°C                                                                                                                                                                                                               |
| <b>Tag</b>               | C-Avi-His                                                                                                                                                                                                           |
| <b>Note</b>              | For research use only                                                                                                                                                                                               |
| <b>Application notes</b> | A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.                                                             |
| <b>Purity</b>            | 98.00%                                                                                                                                                                                                              |
| <b>MW</b>                | 87.80 kDa (predicted); 126.18, 81.32 and 69.37 kDa (reducing conditions)                                                                                                                                            |
| <b>Uniprot ID</b>        | NP_002011.2                                                                                                                                                                                                         |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                                                                      |

### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)