

Product Datasheet

RANKL/TNFSF11/CD254 Protein, Human, Recombinant (rFc) (orb1953388)

Catalog Number	orb1953388
Category	Proteins
Description	RANKL/TNFSF11/CD254 Protein, Human, Recombinant (rFc) is expressed in HEK293 mammalian cells with rFc tag. The predicted molecular weight is 46.18 kDa and the accession number is AAC51762.1.
Tag	N-rFc
UniProt ID	AAC51762.1
MW	46.18 kDa (predicted); 134.21 kDa and 54.28 kDa (reducing condition, due to glycosylation)
Expression System	HEK293 Cells
Biological Origin	Human
Expression Region	A DNA sequence encoding the human TNFSF11 (AAC51762.1) (Gly64-Asp245) was expressed with the Fc region of rabbit IgG at the N-terminus. Predicted N terminal: Ser
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)