

Product Datasheet

PD-L1 Protein, Canine, Recombinant (hFc), Biotinylated (orb1953756)

Description	PD-L1 Protein, Canine, Recombinant (hFc), Biotinylated is expressed in HEK293 mammalian cells with hFc tag. The predicted molecular weight is 51.5 kDa and the accession number is A0A8C0MRV9.
Storage	-20°C
Tag	C-hFc
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	51.5 kDa (predicted); 77 and 145.5 kDa (reducing condition, due to glycosylation)
Uniprot ID	A0A8C0MRV9
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)