

Product Datasheet

SIRP alpha Protein, Human, Recombinant (mFc) (orb1954578)

Description	SIRP alpha Protein, Human, Recombinant (mFc) is expressed in HEK293 mammalian cells with mFc tag. The predicted molecular weight is 63.7 kDa and the accession number is P78324-1.
Storage	-20°C
Tag	C-mFc
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	63.7 kDa (predicted)
Uniprot ID	P78324
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)