

Product Datasheet

CD3E & CD3G Heterodimer Protein, Mouse, Recombinant (hFc & His) (orb1954751)

Description	CD3E & CD3G Heterodimer Protein, Mouse, Recombinant (C-human IgG1-Fc & His & C-human IgG1-Fc) is expressed in HEK293 mammalian cells. The predicted molecular weight is 75.4 kDa and the accession number is NP_031674.1&NP_033980.1.
Storage	-20°C
Tag	C-hFc,His,C-hFc
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	75.4 kDa (predicted)
Uniprot ID	NP_031674.1&NP_033980.1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)