

Product Datasheet

LIF Protein, Human, Recombinant (orb1955622)

Catalog Number	orb1955622
Description	LIF Protein, Human, Recombinant is expressed in HEK293 mammalian cells. The predicted molecular weight is 19.7 kDa and the accession number is P15018-1.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	SDS-PAGE: 99.5%; SEC-HPLC: 98.8%
MW	19.7 kDa (predicted); 23-43 kDa (reducing condition, due to glycosylation)
Uniprot ID	P15018
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)