

Product Datasheet

CCL20 Protein, Human, Recombinant (His) (orb1955788)

Description	CCL20 Protein, Human, Recombinant (His) is expressed in E. coli expression system with His tag. The predicted molecular weight is 10.2 kDa and the accession number is P78556-1.
Storage	-20°C
Tag	N-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	SDS-PAGE: > 90%
MW	10.2 kDa (predicted); 10 kDa (reducing conditions)
Uniprot ID	P78556
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)