

## Product Datasheet

### LIF Protein, Human, Recombinant (His) (orb1956179)

|                          |                                                                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>    | orb1956179                                                                                                                                                                |
| <b>Description</b>       | LIF Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 21.2 kDa and the accession number is P15018. |
| <b>Storage</b>           | -20°C                                                                                                                                                                     |
| <b>Tag</b>               | C-His                                                                                                                                                                     |
| <b>Note</b>              | For research use only                                                                                                                                                     |
| <b>Application notes</b> | A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.                   |
| <b>Purity</b>            | 97.90%                                                                                                                                                                    |
| <b>MW</b>                | 21.2 kDa (predicted); 35-42 kDa (reducing conditions)                                                                                                                     |
| <b>Uniprot ID</b>        | <b>P15018</b>                                                                                                                                                             |
| <b>Expiration Date</b>   | 6 months from date of receipt.                                                                                                                                            |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)