

Product Datasheet

Lymphotactin/XCL1 Protein, Mouse, Recombinant (His) (orb1957077)

Catalog Number	orb1957077
Category	Proteins
Description	Lymphotactin/XCL1 Protein, Mouse, Recombinant (His) is expressed in Baculovirus insect cells with His tag. The predicted molecular weight is 11.5 kDa and the accession number is P47993.
Tag	C-His
UniProt ID	P47993
MW	11.5 kDa (predicted); 20 kDa (reducing conditions)
Expression System	Baculovirus Insect Cells
Biological Origin	Mouse
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)