

Product Datasheet

CCL14 Protein, Human, Recombinant (His) (orb1957412)

Catalog Number	orb1957412
Description	CCL14 Protein, Human, Recombinant (His) is expressed in E. coli expression system with His tag. The predicted molecular weight is 10.2 kDa and the accession number is Q16627-1.
Storage	-20°C
Tag	N-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	10.2 kDa (predicted); 12 kDa (reducing conditions)
Uniprot ID	Q16627
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)