

Product Datasheet

FSTL1 Protein, Human, Recombinant (His) (orb1958514)

Catalog Number	orb1958514
Category	Proteins
Description	FSTL1 Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 34.2 kDa and the accession number is Q12841-1.
Tag	C-His
UniProt ID	Q12841
MW	34.2 kDa (predicted); 47 kDa (reducing condition, due to glycosylation)
Expression System	HEK293 Cells
Biological Origin	Human
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)