

Product Datasheet

LIFR Protein, Human, Recombinant (His) (orb1958620)

Catalog Number	orb1958620
Description	LIFR Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 91 kDa and the accession number is P42702-1.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	SDS-PAGE: 95%; SEC-HPLC: 90.7%
MW	91 kDa (predicted); 125-135 kDa (reducing condition, due to glycosylation)
Uniprot ID	P42702
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)