

Product Datasheet

IGFBP-3 Protein, Human, Recombinant (His) (orb1958757)

Catalog Number	orb1958757
Category	Proteins
Description	IGFBP-3 Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 31 kDa and the accession number is P17936-1.
Tag	N-His
Reconstitution	Reconstituted with sterile deionized water to 0.25 mg/mL.
Purity	96.80%
UniProt ID	P17936
MW	31 kDa (predicted); 40-45 kDa (reducing condition, due to glycosylation)
Expression System	HEK293 Cells
Biological Origin	Human
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)