

Product Datasheet

CEACAM6 Protein, Cynomolgus, Recombinant (aa 35-320, His) (orb1959632)

Description	Carcinoembryonic antigen-related cell adhesion molecule 6 (CEACAM6) belongs to the human carcino-embryonic antigen (CEA) family. Numerous lines of studies have indicated that altered expression of CEACAM6 may have a role in carcinogenesis and development. CEACAM6 Protein, Cynomolgus, Recombinant (aa 35-320, His) is expressed in HEK293 mammalian cells with C-His tag. The predicted molecular weight is 32.73 kDa and the accession number is XP_014979566.2.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Purity	98.00%
MW	32.73 kDa (predicted). Due to glycosylation, the protein migrates to 55-70 kDa based on Tris-Bis PAGE result.
Uniprot ID	XP_014979566.2
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)