

Product Datasheet

ACE2/ACEH Protein, Cynomolgus, Recombinant (His & Avi), Biotinylated (orb1960346)

Catalog Number	orb1960346
Description	ACE2 (Angiotensin I Converting Enzyme 2) is a Protein Coding gene. Diseases associated with ACE2 include Severe Acute Respiratory Syndrome and Neurogenic Hypertension. The protein encoded by this gene belongs to the angiotensin-converting enzyme family of dipeptidyl carboxydipeptidases and has considerable homology to human angiotensin 1 converting enzyme. This secreted protein catalyzes the cleavage of angiotensin I into angiotensin 1-9, and angiotensin II into the vasodilator angiotensin 1-7. ACE2/ACEH Protein, Cynomolgus, Recombinant (His & Avi), Biotinylated is expressed in HEK293 mammalian cells with C-His-Avi tag. The predicted molecular weight is 86.5 kDa and the accession number is A0A2K5X283.
Storage	-20°C
Tag	C-His-Avi
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Purity	98.00%
MW	86.5 kDa (predicted). Due to glycosylation, the protein migrates to 110-140 kDa based on Tris-Bis PAGE result.
Uniprot ID	A0A2K5X283
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)