

Product Datasheet

Latent TGF beta 3 Protein (Primary Amine Labeling), Human, Recombinant (His), Biotinylated (orb1960625)

Catalog Number	orb1960625
Description	Transforming growth factor-betas (TGF- β s) are multifunctional cytokines that have been implicated in the regulation of a broad range of biological processes, including cell proliferation, cell survival, and cell differentiation. And transforming growth factor beta3 (TGFbeta3) is a key protein involved in scar-free healing observed in embryos. Latent TGF beta 3 Protein (Primary Amine Labeling), Human, Recombinant (His), Biotinylated is expressed in HEK293 mammalian cells with N-His tag. The predicted molecular weight is 45.9 kDa and the accession number is P10600-1.
Storage	-20°C
Tag	N-His
Note	For research use only
Purity	98.00%
MW	45.9 kDa (predicted). Due to enzyme lysis and glycosylation, the protein migrates to 43-48 kDa and 50-60 kDa based on Tris-Bis PAGE result.
Uniprot ID	P10600
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)