

Product Datasheet

HDGFRP3 Protein, Human, Recombinant (His) (orb1960932)

Description	Hepatoma-Derived Growth Factor-Related Protein 3 (HDGFRP3) belongs to the HDGF family. HDGFRP3 can be found in testis, heart, spinal cord and brain. HDGFRP3 localizes to the nucleus and contains one PWWP domain. HDGFRP3 enhances DNA synthesis and may have a role in cell proliferation.
Storage	-20°C
Tag	C-6xHis
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	Met1-Thr203
Purity	98.00%
MW	30 KDa (reducing condition)
Uniprot ID	Q9Y3E1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)