

Product Datasheet

IL-13 Protein, Mouse, Recombinant (aa 22-131, His) (orb1961423)

Catalog Number	orb1961423
Category	Proteins
Description	IL-13 Protein, Mouse, Recombinant (aa 22-131, His) is expressed in HEK293 mammalian cells with C-6xHis tag. The predicted molecular weight is 14-30 KDa and the accession number is P20109.
Tag	C-6xHis
Reconstitution	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Purity	100.00%
Protein Sequence	Pro22-Phe131
UniProt ID	P20109
MW	14-30 KDa (reducing condition)
Expression System	HEK293 Cells
Biological Origin	Mouse
Storage	-20°C
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)