

Product Datasheet

TRAIL R2/DR5/TNFRSF10B Protein, Mouse, Recombinant (hFc) (orb1961993)

Description	TRAIL R2/DR5/TNFRSF10B Protein, Mouse, Recombinant (hFc) is expressed in HEK293 mammalian cells with C-Fc tag. The predicted molecular weight is 50-75 KDa and the accession number is Q9QZM4.
Storage	-20°C
Tag	C-hFc
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	Asn53-Ser177
Purity	98.00%
MW	50-75 KDa (reducing condition)
Uniprot ID	Q9QZM4
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)