

## Product Datasheet

### Chicken anti Staphylococcal Protein A, conjugated with Agarose (orb1972344)

<b>Description</b>	Chicken anti Staphylococcal Protein A, conjugated with Agarose
<b>Species/Host</b>	Gallus
<b>Conjugation</b>	Agarose
<b>Tested Applications</b>	IAC
<b>Form/Appearance</b>	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Antibody Type</b>	Primary Antibody
<b>Source</b>	Chickens were immunized with recombinant protein A . Antibodies were isolated from immune eggs, affinity purified on a protein A column and conjugated to N-hydroxysuccinimide (NHS)-activated agarose beads. Beads were washed extensively with after the blocking of the residual NHS sites.
<b>Hazard Information</b>	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
<b>Expiration Date</b>	12 months from date of receipt.

---

#### **Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### **Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)