

Product Datasheet

Chicken anti Rabbit IgG (H+L), absorbed with human and mouse IgG (orb1972356)

Catalog Number	orb1972356
Description	Chicken anti Rabbit IgG (H+L), absorbed with human and mouse IgG
Species/Host	Gallus
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Form/Appearance	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Concentration	See vial for concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Antibody Type	Secondary Antibody
Source	Chickens were immunized with highly purified rabbit IgG isolated from normal pooled rabbit serum, in Freund's adjuvant. Eggs were collected, the IgY fraction prepared and the anti rabbit IgG immunoaffinity purified from a solid phase immunoabsorbant. Anti Rabbit IgG was further absorbed against mouse and human IgG.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)