

Product Datasheet

Chicken Sheep IgG (H+L), conjugated with FITC Antibody (orb1972363)

Catalog Number	orb1972363
Category	Antibodies
Description	Chicken Sheep IgG (H+L), conjugated with FITC Antibody
Species/Host	Gallus
Conjugation	FITC
Reactivity	Sheep
Form/Appearance	Phosphate Buffered Saline, 50% Glycerol with 0.075% sodium azide
Concentration	See vial for concentration.
Purification	Egg-yolk derived IgY Affinity Purified
Tested applications	IHC
Antibody Type	Secondary Antibody
Source	Chickens were immunized with highly purified sheep IgG isolated from pooled normal sheep serum. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were immunoaffinity purified from a solid phase immunoabsorbant and conjugated to fluorescein isothiocyanate.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558