

Product Datasheet

Chicken Soy Protein Antibody (orb1972364)

Catalog Number	orb1972364
Category	Antibodies
Description	Chicken Soy Protein Antibody
Species/Host	Gallus
Conjugation	Unconjugated
Reactivity	Plant
Form/Appearance	Phosphate buffered saline, pH 7.2 with 0.075% sodium azide
Concentration	See vial for concentration.
Purification	Egg-yolk derived IgY Affinity Purified
Tested applications	ELISA, WB
Antibody Type	Primary Antibody
Source	Chickens were immunized with soy (Glycine max) flour protein extract. After multiple immunizations in Freund's adjuvant, eggs were collected and the immunoglobulin fraction prepared. Antibodies were immunoaffinity purified using the soy proteins immobilized on a solid support.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)