

Product Datasheet

ME1 Mouse Monoclonal Antibody (orb1973116)

Catalog Number	orb1973116
Category	Antibodies
Description	ME1 Mouse Monoclonal Antibody
Target	ME1 / malic enzyme 1
Clonality	Monoclonal
Species/Host	Mouse
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Buffer/Preservatives	Phosphate-buffered solution, pH 7.2, 20% Glycerol (pH7). 0.025% ProClin 300 was added as a preservative
Tested applications	IF, IHC, WB
Dilution range	Western Blot 1:500 - 1:1000, Immunofluorescence 1:100 - 1:200, Immunohistochemistry (Paraffin) 1:50 - 1:100
Application notes	Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.
Antibody Type	Primary Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Immunofluorescence of cell. Sample: HeLa, Green: ME1 / malic enzyme 1, monoclonal antibody (orb1973116): 1-100.



orb1973116 ME1 / malic enzyme 1, monoclonal antibody at 1/500 dilution. Lysates at 35 μ g per lane. Nitrocellulose membrane was then blocked for an hour before being incubated with orb1973116 overnight at 4°C.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558