

Product Datasheet

Recombinant mouse CCL19 protein, N-Trx-His (orb1974293)

Catalog Number	orb1974293
Category	Proteins
Description	Recombinant mouse CCL19 protein, N-Trx-His, with a concentration of >2mg/ml, expressed in E.coli, with a molecular weight of 27.3 kDa, and a purity of >90% as determined by SDS-PAGE. It was purified by AC, and has a N-Trx-His tag.
Target	C-C motif chemokine 19
Tag	N-Trx-His
Concentration	>2mg/ml
Purity	>90% as determined by SDS-PAGE
Purification	AC
Protein Sequence	26-108/108
UniProt ID	O70460
MW	27.3 kDa
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Expression System	E.coli
Endotoxins	Not analyzed
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Expiration Date

6 months from date of receipt.



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)