

Product Datasheet

Cosibelimab Antibody (orb1980297)

Catalog Number	orb1980297
Category	Antibodies
Description	Cosibelimab is a fully human anti-PD-L1 monoclonal antibody with a functional Fc domain. It blocks the PD-1/PD-L1 interaction and induces ADCC/CDC, eliminating PD-L1+ cells. It is used in immuno-oncology research for in vitro cytotoxicity assays and in vivo tumor model studies.
Target	PD-1/PD-L1
Clonality	Monoclonal
Conjugation	Unconjugated
Concentration	8.15 mg/mL
Purity	95.00% (May vary between batches)
Antibody Type	Inhibitor Antibody
CAS Number	2216751-26-5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)