

Product Datasheet

Recifercept Antibody (orb1980303)

Catalog Number	orb1980303
Description	Recifercept (TA-46) is a recombinant human soluble fibroblast growth factor receptor 3, a decoy protein that competes for the ligand of mutant FGFR3. Recifercept is used in the study of chondrodysplasia because it restores normal bone growth by preventing negative inhibitory signaling of mutant FGFR3 in patients with pediatric chondrodysplasia. Recifercept is a recombinant human recombinant fibroblast growth factor receptor 3 that competes for the ligand of mutant FGFR3.
Conjugation	Unconjugated
Target	Recifercept
Concentration	7.76 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	98.00%
CAS Number	2226238-09-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)