

Product Datasheet

Pacanalotamab (orb1980321)

Description	Pacanalotamab (AMG 420; BI-836909) is a bispecific T-cell engager (BiTE) that specifically targets BCMA (B cell maturation antigen) and CD3ε. By redirecting T cells to cells expressing BCMA, Pacanalotamab facilitates the T-cell mediated lysis of human multiple myeloma (MM) cell lines [1].
Conjugation	Unconjugated
Target	Pacanalotamab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Purity	98.00%
CAS Number	2251756-52-0
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)