

Product Datasheet

Nimacimab Antibody (orb1980349)

Catalog Number	orb1980349
Category	Antibodies
Description	Nimacimab is a recombinant humanized monoclonal antibody. It can improve diabetes and non-alcoholic steatohepatitis.
Target	Others
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	3.17 mg/mL
Purity	95.40%
MW	144.97 kDa
Antibody Type	Inhibitor Antibody
CAS Number	2098636-09-8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)