

Product Datasheet

Lenvervimab Antibody (orb1980391)

Catalog Number	orb1980391
Category	Antibodies
Description	Lenvervimab (GC1102) is a recombinant human monoclonal antibody of the IgG1 isotype, developed for research into hepatitis B virus (HBV) infection. It is utilized in both in vitro neutralization assays and in vivo preclinical studies to investigate passive immunotherapy against HBV.
Target	HBV
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2055006-79-4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)