

Product Datasheet

Latozinemab Antibody (orb1980394)

Catalog Number	orb1980394
Description	Latozinemab (AL001) is a potent and selective recombinant humanized monoclonal antibody against Sortilin. Latozinemab has a high affinity for Sortilin and interrupts the interaction between the granule protein precursor (PGRN) and the Sortilin receptor. Latozinemab can be used to study of dementia and Alzheimer's disease.
Conjugation	Unconjugated
Target	Latozinemab
Concentration	11.37 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	> 95%
MW	146.28 kDa
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)