

Product Datasheet

Inezetamab Antibody (orb1980412)

Catalog Number	orb1980412
Category	Antibodies
Description	Inezetamab is a bispecific antibody that targets MSLN and agonizes CD40, designed for solid tumor research. It is being investigated in preclinical and clinical studies for its potential to activate immune cells and direct them against tumors expressing mesothelin.
Target	TNF,Mesothelin
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	5.67 mg/mL
Purity	95.00% (May vary between batches)
MW	195.98 kDa
Antibody Type	Inhibitor Antibody
Target Areas	Mesothelin TNF
CAS Number	2568264-76-4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558