

Product Datasheet

Imvotamab (orb1980414)

Description	Imvotamab (IGM-2323), a CD20xCD3 bispecific IGM antibody featuring a dual-action mechanism, promotes physiological T cell activation to avoid both over-stimulation and immune function down-regulation. This compound is under development for addressing B-cell malignant tumors, including multiple myeloma (MM) and non-Hodgkin's lymphoma (NHL) [1] [2] [3] [4].
Conjugation	Unconjugated
Target	Imvotamab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Purity	98.00%
CAS Number	2573121-53-4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)