

Product Datasheet

Flotetuzumab Antibody (orb1980435)

Catalog Number	orb1980435
Category	Antibodies
Description	Flotetuzumab is a CD123/CD3 bispecific DART antibody that activates T cells to eliminate CD123-expressing target cells. It has demonstrated anti-leukemic activity in preclinical models, including AML patient-derived xenografts, supporting its investigation as a potential therapeutic for acute myeloid leukemia in both in vitro and in vivo research.
Target	Immunology/Inflammation related,Integrin
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	2.86 mg/mL
Purity	95.00% (May vary between batches)
Antibody Type	Inhibitor Antibody
CAS Number	1664355-28-5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558