

Product Datasheet

Fidasimtamab Antibody (orb1980440)

Catalog Number	orb1980440
Description	Fidasimtamab (IBI-315; BH2950), a recombinant human IgG1 bispecific antibody, specifically targets, binds to, and inhibits HER2 and PD-1 along with their downstream signaling pathways. Additionally, it bridges PD-1 expressing T cells with HER2 expressing tumor cells, demonstrating potential immunosuppressive and antitumor activities [1].
Conjugation	Unconjugated
Target	Fidasimtamab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	98.00%
CAS Number	2377419-89-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)