

Product Datasheet

Tibulizumab Antibody (orb1980463)

Catalog Number	orb1980463
Category	Antibodies
Description	Tibulizumab (LY 3090106), a tetravalent bispecific monoclonal antibody, selectively binds to B-cell activating factor (BAFF) and IL-17A with dissociation constants (Kd) of 60 pM and 14 pM, respectively. This specificity renders it a valuable tool for autoimmune disease research.
Target	TNF, Interleukin
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	1849636-24-3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)