

Product Datasheet

Valanafusp alfa Antibody (orb1980480)

Catalog Number	orb1980480
Category	Antibodies
Description	Valanafusp alfa (AGT-181) is a recombinant fusion protein combining a chimeric monoclonal antibody targeting the human insulin receptor (HIR) with human iduronidase (IDUA), facilitating brain penetration, and is utilized in researching Mucopolysaccharidosis type I (MPS I).
Target	IGF-1R
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	1815583-32-4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)