

Product Datasheet

Tepoditamab Antibody (orb1980564)

Catalog Number	orb1980564
Category	Antibodies
Description	Tepoditamab (MCLA-117) is a CLEC12A x CD3 bispecific antibody that redirects T cells to eliminate acute myeloid leukemia (AML) stem and progenitor cells. Its application in AML research includes studies of T-cell mediated cytotoxicity in both in vitro and in vivo experimental models.
Target	Integrin
Clonality	Recombinant
Conjugation	Unconjugated
Antibody Type	Inhibitor Antibody
CAS Number	2044679-53-8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558