

Product Datasheet

Zanidatamab Antibody (orb1980569)

Catalog Number	orb1980569
Category	Antibodies
Description	Zanidatamab (ZW25) is a humanized bispecific antibody that simultaneously binds two non-overlapping epitopes (ECD2 and ECD4) of the HER2 receptor. It has demonstrated anti-tumor activity in preclinical research, supporting its investigation in both in vitro and in vivo studies for oncology and HER2-related signaling pathways.
Target	HER
Clonality	Recombinant
Conjugation	Unconjugated
Concentration	3.86 mg/mL
Purity	95.00% (May vary between batches)
Antibody Type	Inhibitor Antibody
CAS Number	2169946-15-8
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558