

Product Datasheet

Andecaliximab Antibody (orb1980587)

Catalog Number	orb1980587
Description	Andecaliximab is a recombinant IgG4 monoclonal antibody targeting matrix metalloproteinase 9 (MMP9). Andecaliximab showed antifibrotic effects in a mouse model of idiopathic pulmonary fibrosis. Andecaliximab can be used to study gastric adenocarcinoma and idiopathic pulmonary fibrosis (IPF), as well as dengue virus.
Conjugation	Unconjugated
Target	Andecaliximab
Concentration	1.93 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	> 95%
MW	150 kDa
CAS Number	1518996-49-0
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)