

Product Datasheet

Bimekizumab Antibody (orb1980645)

Catalog Number	orb1980645
Description	Bimekizumab (Anti-Human IL17A/IL-17F Recombinant Antibody), a humanized monoclonal antibody, selectively neutralizes IL-17A and IL-17F, which promote osteogenic differentiation in human periosteum-derived cells (hPDC), and consequently inhibits inflammation-driven osteogenic differentiation [1].
Conjugation	Unconjugated
Target	Bimekizumab
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	98.00%
CAS Number	1418205-77-2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)