

Product Datasheet

Inotuzumab Antibody (orb1980653)

Catalog Number	orb1980653
Description	Inotuzumab (Humanized Anti-CD22 Recombinant Antibody) is a humanized IgG4 antibody targeting CD22. This naked antibody can be used to synthesize antibody-active molecule conjugators [Inotuzumab ozogamicin] and to study acute lymphoblastic leukemia and non-Hodgkin's lymphoma.
Conjugation	Unconjugated
Target	Inotuzumab
Concentration	6.7 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	> 95%
MW	146.48 kDa
CAS Number	1660159-36-3
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)