

Product Datasheet

Itolizumab Antibody (orb1980662)

Catalog Number	orb1980662
Description	Itolizumab (Anti-Human CD6 Recombinant Antibody) is a recombinant anti-CD6 monoclonal antibody that targets the extracellular SRCR distal domain 1 of CD6, thereby slowing T cell proliferation and inhibiting pro-inflammatory cytokine production. Itolizumab can be used to study conditions such as psoriasis, rheumatoid arthritis (RA), and COVID-19 disease.
Conjugation	Unconjugated
Target	Itolizumab
Concentration	12.76 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Recombinant Antibody
Purity	> 95%
MW	144.82 kDa
CAS Number	1116433-11-4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)